

The Journal of Laryngology and Otology

(Founded in 1887 by MORELL MACKENZIE and NORRIS WOLFENDEN)

EDITED BY

WALTER HOWARTH

WITH THE ASSISTANCE OF

G. H. BATEMAN R. G. MACBETH V. E. NEGUS
R. SCOTT STEVENSON J. F. SIMPSON


1. Original Articles are accepted on the condition that they have not been published elsewhere.
2. Manuscripts should be typewritten, on one side only of the paper, and well spaced.
3. Galley proofs and engraver's proofs of illustrations are sent to the author. Corrections, which should be kept to a minimum, must be clearly marked, and no additional matter added.
4. Blocks illustrating articles are normally paid for by the *Journal*. Authors may, however, be required at the Editor's discretion, to contribute towards the cost of these in cases where this is particularly heavy.
5. Orders for reprints should be sent when returning galley proofs, and for this purpose special forms are supplied.
6. Authors of original communications on Oto-Laryngology in other journals are invited to send a copy, or two reprints, to the *Journal of Laryngology*. If they are willing, at the same time, to submit their own abstract (in English, French, Italian, or German) it will be welcomed.
7. Editorial communications may be addressed to THE EDITOR, *Journal of Laryngology*, c/o HEADLEY BROTHERS, 109 Kingsway, London, W.C.2.
8. The annual subscription commencing with the January, 1950 issue is three guineas sterling, post free, and is payable in advance.
9. Single copies will be on sale at 7s. 6d. each; copies of parts up to Vol. LXIII may be purchased at 4s. each.
10. All subscriptions, advertising and business communications should be sent to the publishers, HEADLEY BROTHERS, 109 Kingsway, London, W.C.2.

United States of America

Orders for this *Journal* may be sent either through a local bookseller, or to Messrs. G. E. STECHERT & Co., 31-33 East 10th Street, New York, or they may be sent direct to the publishers, HEADLEY BROTHERS, 109 Kingsway, London, W.C.2, England.

Please mention *The Journal of Laryngology* when replying to advertisements

THE AMPLIVOX AUDIOMETER



The Basic Instrument of Modern Otology

NEW TESTS NOW POSSIBLE

Speech-loss measurement in decibels now obtainable through dynamic microphone and speech input meter.

Loudness balance control gives measurement of the recruitment factor.

Input for group-test equipment and gramophone.

CONSISTENT ACCURACY OF RESULTS

Eleven exact test frequencies 125-12,000 c.p.s.

Voltage compensation ensures that signal intensity does not vary should mains voltage do so.

CALIBRATED MASKING TONE

Masking intensity varies equally with normal hearing acuity over range 125-4,000 c.p.s., permitting accurate masking control. Calibrated in decibels.

SIMPLICITY OF OPERATION

Twin matched headphones (on double headband) for all air conduction tests; bone conductor on single headband, or interchangeable with one 'phone.


Hearing Loss dial calibrated in 5 dB steps from minus 10 dB to plus 100 dB.

Same figures on dial are read for air and bone conduction, and for all frequencies.

Full details are available from the manufacturer, who will gladly arrange demonstrations, if required.

AMPLIVOX LTD., 2 BENTINCK ST., LONDON, W.I. (Welbeck 2591)


Please mention *The Journal of Laryngology* when replying to advertisements





WATER

SUCTION APPARATUS

Particulars on application



FOOT POWER



ELECTRIC

DOWN BROS. & MAYER & PHELPS Ltd.

Correspondence: 23 PARK HILL RISE, CROYDON

Telephone: Croydon 6133, 6134, 6135

Showrooms: 32/34 NEW CAVENDISH STREET, LONDON, W.1


Telephone: WELBECK 3764

DOWN BROS. AND MAYER & PHELPS HAVE AMALGAMATED.
THE PERSONAL MANagements REMAIN AS HERETOFORE.

Please mention *The Journal of Laryngology* when replying to advertisements

PRESCRIBE THE EASY EFFECTIVE METHOD

When spraying, the patient can hold his head in a natural position, yet the solution spreads throughout the nasal cavity, assuring thorough treatment. Effective dropper-application, on the other hand, requires a complicated posture technique which patients are apt to neglect in home treatment.



DROPPER ?



**DeVilbiss
ATOMIZERS**



ATOMIZER ?

These illustrations, which are based on X-ray research, clearly demonstrate the advantages of the atomizer as compared with the dropper
 The **AEROGRAPH Co., Ltd., Lower Sydenham, London, S.E.26**
 Sole Distributors of DeVilbiss Products in the United Kingdom.
 Telephones: Sydenham 6060—8 Lines.

CVS—225

For Personal Protection

VAPEX

The pleasant, unobtrusive aroma of Vapex, its efficacy in the relief of the common cold and as a prophylactic, make it particularly suitable for personal use by the medical profession.

<p>INHALANT</p> <p>A drop on the handkerchief or on the end of the pillow.</p>	<p>INHALER</p> <p>Convenient and unobtrusive for use during the day. May be easily recharged.</p>	<p>VOLATOL</p> <p>For topical application in chest colds. Non-greasy and does not stain.</p>
---	--	---

THOMAS KERFOOT & Co. Ltd., Vale of Bardsley, Lancashire.

V260

Please mention *The Journal of Laryngology* when replying to advertisements



THE S.D. PORTABLE AUDIOMETER A.C./D.C.

THIS genuinely portable Audiometer, built into a strong carrying case, offers great accuracy combined with convenience in use.

It incorporates an extremely stable beat frequency oscillator giving a continuous range of frequency, and a built-in valve voltmeter not only allows accurate zero setting of frequency but ensures a constant result whatever the conditions of voltage or of valve life.

A crystal microphone and amplifier give excellent speech reproduction and in addition enable speech level tests to be made. A masking note of mixed low frequency and a manual mute are also provided.

Our Peepshow for testing the hearing of young children can be used with this Audiometer.

ALFRED PETERS & SONS LTD.
FARGATE HOUSE · SHEFFIELD 1

Please mention *The Journal of Laryngology* when replying to advertisements



THE SD/2 AUDIOMETER

THIS instrument offers not only great accuracy but incorporates some unique features, which allow a most comprehensive range of tests to be made and ensures a constant result under all conditions of voltage and valve life.

It includes among many facilities an automatic interruptor circuit, as well as a continuous tone: twin matched earphones allowing the test and masking tones to be switched to either ear, while a balance control enables the response of each ear to be matched simultaneously.

A built-in crystal microphone and amplifier give excellent speech reproduction and allow speech level tests to be made, while the built-in meter and cathode ray tube ensure continuing accuracy of frequency and output.

Our Peepshow for testing young deaf children may be used with this Audiometer.

ALFRED PETERS & SONS LTD.
FARGATE HOUSE · SHEFFIELD 1

Please mention *The Journal of Laryngology* when replying to advertisements

CONTENTS

RESIDUAL CAPACITY TO HEAR OF PUPILS IN SCHOOLS FOR THE DEAF. A. I. Goodman (Manchester)	551
REMARKS ON A THEORY OF THE EVOLUTION OF THE TETRAPOD MIDDLE EAR. F. R. Parrington (Cambridge)	580
MID-LINE CONGENITAL MALFORMATIONS OF THE NOSE. A. S. Handousa Bey (Cairo, Egypt)	596
PHARYNGO-ŒSOPHAGEAL DIVERTICULUM. S. Horowitz (Glasgow)	600
CLINICAL RECORD—	
Two Cases of Lateral Sinus Thrombosis. E. H. Hadfield and R. T. Hinde (Oxford)	605
CLINICAL NOTE—	
Individual Inserts for Hearing Aids. E. A. Dennison, A. C. Deverell, T. S. Littler and E. Whetnall (London)	610
SOCIETIES' PROCEEDINGS—	
Royal Society of Medicine—Section of Laryngology	614
ABSTRACTS—	
Ear	629
The Relation between Hearing Loss for Specific Frequencies and the Distance at which Speech can be Identified. E. Thayer Curry; Chorda Tympani Nerve Section and Tympanic Plexectomy. Samuel Rosen.	
Larynx	629
Chondroma and Chondrosarcoma of the Larynx. Melvin R. Link; Rehabilitation of the Larynx in Cases of Bilateral Abductor Paralysis. De Graaf Woodman.	
Œsophagus	630
Œsophageal Varices. Charles B. Ripstein.	
Miscellaneous	631
Vertebrate Magnets for Removal of Foreign Bodies from the Air and Food Passages. Chevalier Jackson and Chevalier L. Jackson; Tantalum in Rhinoplastic Surgery. Samuel L. Fox; Mitral Stenosis Associated with Left Recurrent Nerve Paralysis. D. F. Lewis and J. W. Scott; Anticoagulant Therapy. C. H. Jaimet; Postgraduate Training in Oto-laryngology. Francis L. Lederer and Joseph G. Schoolman; Functional and Anatomic Relation of Sphenopalatine Ganglion to the Autonomic Nervous System. David Higbee.	
GENERAL NOTES	634

For Advertisement space in this Journal apply to :
HEADLEY BROTHERS, 109 Kingsway, London, W.C.2

Please mention *The Journal of Laryngology* when replying to advertisements