Standard Setting s13

This process repeated in the third round. Statements attaining consensus moved to the final report. The remaining statements did not attain consensus.

Results/Outcomes: Seventy-seven experts participated in the first, 75 in the second, and 74 in the third round. Twenty-three statements attained consensus. Twenty-one statements did not attain consensus.

Conclusion: A Delphi technique was used to establish consensus regarding the numerous complex factors influencing treatment of the E&ME in RSEs. Twenty-three statements attained consensus and can be incorporated into guidelines to advance the EMT treating the E&ME.

Prehosp. Disaster Med. 2025;40(Suppl. S1):s12-s13 doi:10.1017/S1049023X25000494

The Use of Point-of Care Diagnostic Techniques in an EMT1 Fixed Improves Diagnostic Accuracy and the Detection of Seriously III Patients

Alberto Hernández-Tejedor MD, PhD [0], Sara Isabel Montero Hernández RN, David Caballos Villar MD, Víctor Sainz Ruiz de León

SAMUR-Protección Civil, Madrid, Spain

Background/Introduction: The use of diagnostic technology at the bedside is increasingly common in hospitals and is spreading to EMS.

Objectives: To detect the contexts in which the use of point-of-care blood test and ultrasound is feasible and provides benefit in an EMT1 fixed.

Method/Description: Review of the procedures of the SAMUR-Protección Civil and the bibliography that supports them. Summary of the publications carried out by our group. Results/Outcomes: The basic blood analysis with blood gases, biochemistry and electrolytes allows to detect electrolyte disturbances secondary to gastrointestinal infection, diabetic decompensation, and others, to estimate the severity of a trauma patient that is not clearly unstable, to guide resuscitation in shock, to study the patient with chest pain, to assess organic involvement in infectious conditions, etc.

The use of ultrasound at the bedside has shown to be a great diagnostic aid in many pathologies: trauma (frequent in the context of a catastrophe), COVID-19 (detection of infiltrates), cardiorespiratory arrest (reversible causes), pulmonary embolism (high risk in people trapped), study of collections and abscesses, heart disease, etc. and it is useful for multiple techniques such as vascular cannulations (peripheral and central), orotracheal intubation, and drainage.

The ESP EMT1 SAMUR-PC is equipped with ultrasound and point-of care blood analysis, as well as our ALS units in Madrid.

Conclusion: The inclusion of ultrasound and blood tests among the diagnostic capabilities of an EMT1 fixed can

improve the ability to detect serious pathology and guide initial treatment.

Prehosp. Disaster Med. 2025;40(Suppl. S1):s13 doi:10.1017/S1049023X25000500

Evaluating the Effectiveness of Dry Decontamination Methods for Mass Decontamination: A Scoping Review of Current Research

Eman Alshaikh MD, MS, FIBODM^{1,2}, Attila Hertelendy PhD¹, Fadi Issa MD^{1,2}, Amalia Voskanyan RN¹, David Digregorio PA, MSEM, MPAS¹, Gregory Ciottone MD^{1,2}

- Beth Israel Deaconess Medical Center (BIDMC), Boston, MA, United States
- 2. Harvard Medical School, Boston, MA, United States

Background/Introduction: In the event of CBRN mass casualty incidents, a specialized approach is essential that includes decontamination procedures. There are two types of decontamination systems most often used for CBRN incidents: wet and dry. While wet decontamination is more common, it can be problematic in cold weather or when logistical equipment is unavailable. This scoping review seeks to assess the effectiveness of dry decontamination.

Objectives: Evaluating the effectiveness of dry decontamination methods for mass decontamination.

Method/Description: A scoping review was done on dry decontamination using various databases, including MEDLINE, CINAHL, Embase, Web of Science, and Scopus. Following the PRISMA approach, nine eligible articles were included in this scoping review. The methods of dry decontamination, types of decontamination materials, methods of analyzing decontamination, and the main conclusions from each study were extracted, summarized, and compared.

Results/Outcomes: The review found that dry decontamination is an effective method (4 articles), particularly when dealing with liquid contaminants (1 article). However, its efficacy diminishes in cases of more hairy parts of the body (1 article) and when dealing with chemicals in particle form. Effectiveness of dry decontamination varies depending on the material used and the chemical agent involved (3 articles). Moreover, the technique of blotting and rubbing was found to increase effectiveness (1 article).

Conclusion: This scoping review found that dry decontamination can be effective and may be an alternative to wet decontamination in CBRN mass casualty incidents, particularly for liquid contaminants. Its effectiveness varies depending on factors such as the type of contaminants, and the decontaminant material used.

Prehosp. Disaster Med. 2025;40(Suppl. S1):s13 doi:10.1017/S1049023X25000512

